


## DIRECTIONS FOR USE:

**Allow the sealed foil pouch with the word " 1 Test", the bag labeled "Extraction Tube", the bag labeled "SARS-COV-2 Buffer", and the "Sterilized Swab" to warm up to room temperature (60-80 degrees fahrenheit).**



1. Remove the SARS-COV-2 rectangular test device from the sealed foil pouch that says 1 Test and use it as soon as possible.



2. Place the test device on a clean and level surface. Best results will be obtained if the attempt is performed immediately after opening the foil pouch.
3. Unscrew the whole cap of the Extraction Tube Nozzle.



- 4.
5. Take out on (1) bottle from the SARS-COV-2 Buffer,



Remove the bottle cap, and add all the liquid into the Extraction Tube. (See Illustration #1 below.)

6. Open the Swab container, extract the Swab and insert it into your nose. Swirl it around 10x in each nostril.



7. Place the sterilized swab specimen in the Extraction Tube (which already contains the SARS-COV-2 Buffer). Rotate the Swab for approximately 10 seconds while pressing the head against the inside of the Extraction Tube. (This will release the antigen in the Swab). (See Illustration #2 below).
8. Remove the Swab after squeezing the Swab head against the inside of the Extraction Tube (in order to expel as much liquid as possible from the Swab. (See Illustration #3 below). Discard the Swab in accordance with your biohazard waste disposal protocol.
9. Screw on and tighten the cap on the Extraction Tube, and then shake the Extraction Tube vigorously to mix the specimen and the SARS-COV-2 Buffer liquid. (See Illustration #4 below).
10. Add 3 drops of the solution in the Extraction Tube (approximately 80ul) to the sample well, leave it flat on the table and then start the timer. (see Illustration #5 below). Read the result 10-20 minutes later. Don't interpret the result after 20 minutes.

